

ASR Mobile Cart Best Practices

Infection Prevention

Patient and Caregiver Safety

Patient Privacy

Patient Experience

Inspect at Start of Shift

- Charged up
- Computer up
- Scanner working
- Badge Reader working and Log In
- Check Connectivity

Wipe all surfaces down (sani-wipe) worksurface, handles, keyboard, keyboard tray, mouse, mouse pad area, barcode scanner
Wipe down Monitor Exterior (Note: Never wipe monitor screen with Sani Wipe or Alcohol Wipes...use water on Microfiber cloth)

Check Your Supplies

- Common Supplies - Fill what is needed, excess supplies add more weight making Car harder to push
- Special Supplies - Add what is needed
- Instruments needed for the shift

During the Shift

- Plug in to top off battery when cart is not in use - even for 5 minutes extends run time
- Wipe down the cart after leaving each room
- Keep cart out of main traffic areas in hallways
- Watch where you are going with cart - adjust monitor to provide clear view of your intended direction
- Be Careful with cart going over thresholds into carpeted area or returning to hard surface flooring
- Be very careful when entering and exiting elevators - some do not always align correctly with the floor
- Turn off display or go to main screen when leaving Cart for long period of time
- Face and engage patient / Family when Charting to increase Patient Experience

Issues

- Try to resolve - sometimes a simple reboot works, plug the Cart in,
- Ask a co-worker for quick help
- Report immediately to IT (open a Ticket)
- Connectivity or Signal Drops - report to IT



For additional information or help, contact us at sales@asrhealthcare.com or 877.272.8738

"we know Point-of-Care"

ASR Healthcare - 4319 Abbotts Bridge Road Suite 5 Duluth, GA 30097